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
- The Hovercraft is designed to help students experience frictionless motion, thus better understand Newton's Laws.
- Its large platform provides enough area for the rider to comfortably sit while riding.
- The durable Nylon skirt will withstand the rigors of the classroom environment.
- The optional Cordless Air Source is a convenient way supply air to the Hovercraft, however most vacuums/blowers used in a wood shop will also be sufficient to lift the Hovercraft.

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-PASCO SE8747

Key Features

- The Student "Becomes" the Dynamics Cart
- Experience the Feel of Physics Concepts

When a student sits on top of the Kinesthetics Cart, he or she actually becomes part of the experiment.

Active involvement makes for "fun physics," but, more importantly, it provides a sensory experience, or "kinesthesia," that teaches what a concept such as Newton's First Law "feels like."



Students' common misconceptions are quickly eradicated.

The experience is retained in "muscle memory" and helps the student understand and remember the concepts.

一維動力台車 (1D Kinesthetics Cart)

-PASCO SE8747

Features:

1. Ball bearings mounted on the bottom of the upper cart allow it to slide off the lower cart
2. Low "Roller Blade" Wheels on both carts allow for smooth, constant travel
3. Coupler keeps carts piggy-backed until rider removes it
4. Cart length of 90.5 cm

Typical Experiments of 1D Kinesthetics Cart

1. Newton's First Law
2. Newton's First Law
3. Newton's Second Law
4. Newton's Second Law
5. Coin-toss Misconception Demo*
6. Projectile Motion*
7. Uniform *vs.* Accelerated Motion*
8. Newton's Third Law
9. Center of Mass -- The Boardwalk*
10. Simple Harmonic Oscillator*
11. Human Oscilloscope*
12. Collisions: Newton's Third Law
13. Newton's Third Law Misconceptions
14. Newton's Third Law Misconceptions
15. Newton's Third Law Misconceptions


Experimental Manual

*These experiments require additional equipment not included with the Kinesthetics Cart.



Sample: Experiment 1 - Newton's First Law Demonstration


1. One rides on the piggy-backed carts at initial velocity, V_o . He removes the coupler as he approaches the block fastened to the floor.
2. The lower cart strikes the block and stops, but the upper cart and student continue to move forward at velocity, V_o .
3. The slight jolt that the student feels when the wheels hit the floor will help ingrain Newton's First Law in the student's memory.
4. **Newton's First Law:** An object in motion continues in motion, unless acted upon by an external force.

Experiment 6


-Independence of the X and Y Motion of a Projectile

Experiment 9: Center of Mass Conservation of Momentum

Newton's Third Law- "A Crushing Experience"

Experiment 6a: Independence of the X and Y Motion of a Projectile



Fig. 6-1: A ball is launched vertically from a stationary cart (stationary observer).

Experiment 4: Newton's Second Law -“The Bucket Accelerator”


二維動力台車(2D Kinesthetics Cart)

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- With this circular, 6-wheeled cart, students can "feel" the forces and accelerations that occur during two-dimensional motion such as circular and trajectory motion. This cart enables students to directly experience centripetal force, and it allows quantitative verification of the relationship between centripetal force, orbital period and the radius. As they sit on the cart, students will finally be convinced that centripetal force is directed radially inward.
- Manual-2D KINESTHETICS CART

The Foucault Pendulum Demonstration

The Foucault
Pendulum in action.

In the rotating frame of the rider, the plane of the pendulum rotates, although in the lab frame its plane of oscillation remains fixed.


雲霄飛車實驗 (Complete Roller Coaster)


The Roller Coaster including the track and Mini Cars and some accessories allows the quantitative study of energy and motion.

Key Features:

1. 3-Car Roller Coaster
2. Quantitative Studies of Energy Conservation
3. Easy to Change Track Configurations

雲霄飛車實驗(Complete Roller Coaster)

Typical Applications:

- 1. Conservation of Energy**-Release the Mini Car and measure its velocity and height at several points along the track. Use these values to calculate total energy of the Mini Car. Frictional losses are less than 5%.
- 2. Constant Acceleration**-Several straight inclined sections can be used to measure and demonstrate constantly accelerated motion.
- 3. Projectile Motion/Conservation of Energy** -- Use the initial height of the Mini Car to determine its speed as it flies off the end of the track. Using this speed and height above the ground when it leaves the track, predict where the Mini Car will land.

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
- **Experiment 1: Projectile Motion**


- The purpose of this experiment is to predict and verify the range of a ball launched at an angle.
- The initial velocity of the ball is determined by shooting it horizontally and measuring the range and the height of the Launcher.

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Double Cone rolls upward on a plane

- **Demonstration: Ascending Cone**-Defy gravity as our cone appears to roll upward.
- This visual demonstrator of the center of gravity consists of a wooden frame with diverging rails and a double-ended cone.
- The cone appears to roll upward; in actuality, it is the center of mass that is moving downward.

The physical explanation:


(a) **The center of mass of the cone is descending.**

(b) The motion may be described by the energy conservation law for the cone-Earth as follows:

$$(1/2)mv^2 + (1/2)I\omega^2 + mgw = \text{constant}$$

$I = (3/5)mR^2$ = rotational inertia of double cone
with respect to symmetry axis

Where m = mass (g), R = radius (cm), v = velocity (m/s²) and ω = angular velocity (rad/s), w = instantaneous height of center of mass of double cone over zero level of potential energy.

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Maxwell Wheel

- [Refer ppt file](#)
- Useful for demonstrating momentum, conservation of energy and torque, this apparatus is basically a large metal flywheel suspended by two strong cords in a metal frame.
- The cords are wound around the shaft of the wheel which is then released. The wheel will unwind as it falls but will wind itself back up as the momentum carries it upward in the opposite direction.
- This oscillation process will continue for several moments before gradually coming to a stop as the wheel slowly loses momentum and travels less each time.

Air Cannon (空氣砲彈) or Airsooka Air Gun

-A great demonstration of the energy that can be stored in waves.

- The Air Cannon uses a **vortex of air** for ammunition.
- Its unique shape creates a **stable toroidal vortex**. Pull back the flexible membrane, release and the invisible wave front of air can hit a target 20 feet away!

Airzooka air gun

- Launch a full air assault without ever leaving the ground with the 'Airzooka air gun'.
- Airzooka is a hand held fun gun that fires a ball of air with enough power to mess up a persons hair, ruffle their shirt or scatter papers from a distance of 20 feet away.
- By pulling back the membrane and releasing - a shock wave of air is sent hurtling at your target.
- Requiring no batteries or electricity, Airzooka works with an elasticated air launcher.
- Pull back and release the launcher like a catapult, firing the ball of air out of the air gun.
- Because it shoots air, you never run out of ammo!


The Vapor Ring structure is called a "Toroidal Vortex".
 Toroidal = doughnut shaped
 Vortex = whirling or circular motion of a liquid or air, tending to form a low pressure area in the center.
 Tomadoes, hurricanes, and whirlpools are naturally occurring vortices.

Vapor-Ring Physics


QUESTION; How does The Vapor-Ring (TOROID) move forward?
 Think of the cross-section of rotating air in the Vapor-Ring like the tire on a tractor moving over the ground. The friction of the rotating air on the outside surface of the Vapor-Ring against the surrounding stationary air makes the ring tractor through the air. There is an overall circular flow of air. Pressure is not created in front of the Vortex. The Vapor in the Vortex makes it appear that air is moving from left to right. The backward flow of air cancels out the forward flow of air.

Air flow is slowed at edge of opening causing a vortex to form with the fast moving air through the center ("Coanda effect")

Low pressure at center of Toroidal Vortex holds Vapor around center of ring

stationary air creates friction against the rotating air in the vapor-ring

Air friction makes the Vapor-Ring grow larger in diameter and pulls Vapor away from the ring as it moves forward.

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
Step 1: Use your thumb and forefinger to support the chain in the center of the ring. At the same time, your middle finger is close to your thumb.

Step 2: Let the ring lean against your middle finger loose the ring, then the ring will be hitched by the chain.

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- 伽利略溫度計(Galileo Thermometer)
- 不停地喝水的鳥(Thermodynamic drinking bird)
- 熱動力沸騰器(Thermodynamic Boiler)-液氣相轉變示範實驗
- Hand Bubbler Boiler Pen Demonstrates Energy Transfer
- 輻射計(Radiometer)-光熱轉輪
- 咖啡杯上的史特林引擎(Stirling engine)
- 翅膀揮舞的蝴蝶-是熱縮冷脹嗎?
- 紫外光偵測珠(Ultraviolet Detecting Beads)
- 太陽能袋World's Largest Solar Bag

Galileo Thermometer (伽利略溫度計)

Be invented by Galileo Galilei four centuries ago
To indicate the approximate current temperature.

How does a Galileo thermometer work?

- Based on a thermoscope invented by Galileo Galilei in the early 1600s.
 - ⇒ A simple, fairly accurate thermometer, today it is mostly used as decoration.
- Be consisted of a sealed glass tube that is filled with water and several floating bubbles.
- ✓ The bubbles are glass spheres filled with a colored liquid mixture. This liquid mixture may contain alcohol, or it might simply be water with food coloring.
- ✓ Attached to each bubble is a little metal tag that indicates a temperature. A number and degree symbol are engraved in the tag.

- **These metal tags are actually calibrated counterweights.**

The weight of each tag is slightly different from the others.

- ✓ Since the bubbles are all hand-blown glass, they aren't exactly the same size and shape. **The bubbles are calibrated by adding a certain amount of fluid to them so that they have the exact same density.**
- ✓ So, after the weighted tags are attached to the bubbles, each differs very slightly in density (the ratio of mass to volume) from the other bubbles, and the density of all of them is very close to the density of the surrounding water.

- An object immersed in a fluid experiences two major forces: **the downward pull of gravity and the upward push of buoyancy.** It is the downward force of gravity that makes this thermometer work.
- The basic idea is that as the temperature of the air outside the thermometer changes, so does the temperature of the water surrounding the bubbles.
- As the temperature of the water changes, it either expands or contracts, thereby changing its density. So, at any given density, some of the bubbles will float and others will sink.
- **The bubble that sinks the most indicates the approximate current temperature.**

Galileo thermometer

- Consider the example: Let's say there are five bubbles in the thermometer:


1. A blue bubble represents 60°F
2. A yellow bubble represents 65°F
3. A green bubble represents 70°F
4. A purple bubble represents 75°F
5. A red bubble represents 80°F

- The blue bubble (60°F) is the heaviest (densest) bubble, and
- each bubble thereafter is slightly lighter, with the red bubble being the lightest.
- Now, let's say the temperature in the room is 70°F.
- Since the surrounding air is 70°F, we know the water inside the thermometer is also about 70°F. The blue and yellow bubbles (60 and 65°F, respectively) are calibrated so that they have higher densities than the water at this temperature, so they sink.
- The purple and red bubbles each have a density that is lower than the surrounding water, so they float at the very top of the thermometer. Since the green bubble is calibrated to represent 70°F, the same temperature as the water, it sinks slightly so that it is floating just below the purple and red bubbles -- thereby indicating the room's temperature!

Solution

- A Galileo thermometer combines [Archimedes' principle](#) with the fact that [liquids generally expand faster](#) with increasing temperature than solids do.
- Each sphere in the thermometer has [an average density \(a mass divided by volume\)](#) that is very close to that of the fluid in the thermometer.
- As stated in Archimedes' principle, if the sphere's average density is less than that of the fluid, the sphere floats and if the sphere's average density is more than that of the fluid, it sinks.
- But the fluid's density changes relatively quickly with temperature, becoming less with each additional degree.
- Thus as the temperature of the thermometer rises, the spheres have more and more trouble floating.
- Each sphere's density is carefully adjusted so that it begins to sink as soon as the thermometer's temperature exceeds a certain value.
- At that value, the expanding fluid's density becomes less than the average density of the sphere and the sphere no longer floats.
- The spheres also expand with increasing temperature, but not as much as the fluid.
- Here is a picture of a combined Galileo thermometer and [simple barometer](#).
- In addition to measuring the temperature with floating spheres, this device measures the outside air pressure with a column of dark liquid.
- It has a trapped volume of air that pushes the liquid (visible at the bottom of the unit) up a vertical pipe when the outside air pressure drops.

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- 一個熱能與機械能互相轉換的熱力學過程。
- The dunking is due to the evaporating water, which cools the head reducing the pressure of the gas inside. The liquid inside the bird rises up from his tail towards his head, causing the bird to tip and appear to be drinking. As the bird dunks, the liquid returns to its tail, and the process begins again.

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

Thermodynamic Drinking Bird

- The 'Drinking Bird' is disguised as a great demonstration of thermodynamics.
- When the head of the bird is moistened with water it will repeatedly tip over and appear to be drinking from a cup placed in front of it.
 - The dunking is due to the evaporating water, which cools the head reducing the pressure of the gas inside. The liquid inside the bird rises up from his tail towards his head, causing the bird to tip and appear to be drinking. As the bird dunks, the liquid returns to its tail, and the process begins again.
 - In essence this is a small, cleverly designed thermodynamic engine which only needs water to run!
 - Avoid breaking glass, contains methylene chloride (二氯甲烷) and may stain clothing and other surfaces.
 - ✓ Do not swallow or allow to come into contact with eyes. Flammable. Keep away from flame and other heat sources. CHOKING HAZARD:

- **Setup:** Wet cloth which covers the tube in the bird (i.e. wet bird's head) and set up the bird with a beaker in front of it. The gas in the birds, which consists of various organic fluids, will condense causing the birds to "drink". **The water cools and reverses the condensation process** and the bird settles back down.
- **Concepts Displayed: condensation, expansion, Rankine cycle**
- [Click Here to see a short videoclip of the demonstration.](#)
- [Click here for a more detailed explanation of how it works](#)

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
Hand Bubbler Boiler Pen

Demonstrates Energy Transfer

- Hold the bottom chamber between your fingers and watch as your body heat warms the liquid.
- As it starts to vaporize, it expands moving it's way through the tubes to the upper chamber. When all the liquid has reached the top, a bubbling effect is created!
- Let go and the liquids cool, and move back down to the lower chamber. Great for science class! Includes matching stand.
- **WARNING:** Not intended as a toy for children.
 1. **Contains Ethyl Alcohol.** Avoid contact with eyes. In case of eye contact, flush thoroughly with water. If irrigation persists, get medical attention.
 2. Do not use near heat or flames. Only warm using your hand. Liquid will stain if glass is broken.
 3. Keep out of reach of children.

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
Ultraviolet Detecting Beads - UV Science and Fun

- These UV beads are an inexpensive fun way to make the invisible become visible.
- The solar pony beads are pale in visible light, such as indoors or if they are shielded by a UV screen such as automobile glass or sunglasses. The Ultra-violet sensitive beads contain a pigment which changes color when exposed to UV light from the Sun or other UV source.
- Many of today's space telescopes observe colors of light not in the visible spectrum such as infrared and ultraviolet.
- One can make a simple bracelet from a standard pipe cleaner or strip of rawhide which can function as a UV detector. Use your imagination to think of other things to use them for.
- Many other experiments are possible as well... [test sunglasses for the effectiveness of their UV protection](#), or even [suntan lotion](#).
- In five assorted colors, all beads are pale, off-white in the absence of UV light.
- The ultraviolet beads will cycle back and forth (to bright colors and back again) over 50,000 times!
- Become an ultraviolet detective with these amazing UV-sensitive beads! It's solar magic in action!
- Make bracelets, test the effectiveness of sunscreen, or experiment with color while you learn about the power of sunlight. Included is a science activity guide with lots of fun experiments with solar radiation and light.

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Ultraviolet Detecting Beads - UV Science and Fun

- Become an ultraviolet detective with these amazing UV-sensitive beads! It's solar magic in action! Make bracelets, test the effectiveness of sunscreen, or experiment with color while you learn about the power of sunlight. Included is a science activity guide with lots of fun experiments with solar radiation and light.
- Ultraviolet Beads glow with bright colors when exposed to the sun or other UV light source, but remain off-white when indoors or not in the presence of UV.
- Students can use UV Beads to design and conduct experiments that test the effectiveness of UV protection: sunglasses, sun-block, and automobile glass for example.
- Combine science and art to create jewelry and other fun art projects using the beads. String the colorful little UV indicators together to create necklaces or bracelets that double as UV detectors. " "

World's Largest Solar Bag

- This amazing item demonstrate principles of density and buoyancy, as well as solar energy and laws of gases.
- The largest on the market today, this solar bag is 50 feet long and 3 feet in diameter and can be cut in half to make two bags the same size or two different sizes.
- The thin, black plastic bag holds over 350 cubic feet of air.

Experiments (實驗步驟) :

- Fill the bag by running, then tie up the ends and let the solar energy of the sun heat the air inside the bag.
- Tether the bag with string and watch it rise into the sky in only five minutes.

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- 韋氏起電機(Wimshurst's Electrostatic Generators)
- 范式(范德格拉夫)起電機(van de Graaf Generator)
- 靜電展示與應用
- 比比哪個磁力比較強？
- Spinning Magnet, Levitating Vortex 物體騰空懸浮歷久不墜
- Levitron-Floating Gyroscope

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Wimshurst's Electrostatic Generators

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


Demonstration:


- Hair Raising
- Electric Wind
- Lightning
- St. Elmo's Fire


Schematic shows positive charges carried upward from Teflon pulley to collector sphere

Start-Up

Combs reach ionization intensity


Dome reaches ionization intensity

Demonstrations: Hair Raising

<http://www.sciencefirst.com/>


Instructions file:

To raise hair, stand on footstool or other plastic (nonconducting) platform to insulate yourself.

Demonstration: Electric Wind

- **By attaching the conductor or needle, you have created an electric wind.** The conductor is in the form of a sturdy, light, thin metallic rod six to eight inches long (for instance, a darning needle) on the body of the collector dome, radially outwards. Use tape or clay to attach.
- The concentration of charges at the tip of the needle will be so intense that it will ionize air in its neighborhood. Negative ions will rush towards the collector dome and neutralize their charges. Positive ions, however, move away (due to electrostatic repulsion) from the generator and do not get neutralized.
- As the generator is continuously running, it keeps supplying more and more positive ions at a fast speed. The ions running away form a wind called “electric wind” which blows away (radially outward) from the generator.

Using tape or clay to attach a conductor or needle :

- **Generate statics**
- **Turn a vane**
- **Spin a spinner**
- **Deflect a flame**
- **Rotate the collector dome**

Demonstration: Electric Wind

- **Generate Statics:** The wind is strong enough for its effects to be experienced as far away as 0 feet from the generator. It may not deflect a flame that far away but will certainly impart statics to your clothing which would cling to your body; or to a paper that would cling to your hand or to the wall.
- **Turn a vane:** Place a vane, such as a child's pinwheel, in front of the conductor.
- It will turn in the direction of the wind. See for yourself what the wind direction is and see if you can form some idea of how strong the wind is. Try a vane that is slightly stiff and requires a stronger wind to turn it.
- **Spin a spinner:** Make a small spinner using aluminum foil " across with 4 - 6 blades. Use a sharp pin to act as axis for spinner and mount the pin on a wooden or plastic stick. Try placing 2 beads on each side to localize the spinner. When brought near the conductor, the electric wind will spin it.
- **Deflect a flame:** Bring a lighted candle near the conductor. The flame is deflected away from the generator in the same manner as an air draft.

Demonstration: (2) Electric Wind

- **Rotate the collector dome:** Show how an actual (electric) wind can be created by ionized air molecules running away from the pointed conductor. The ionized molecules move away from the sharp or rounded end of the conductor in great numbers and at great speeds. This, according to **Bernoulli's Principle**, produces a low-pressure region in front of the tip of the conductor. The rear end of the conductor (attached to the dome) remains at normal pressure. This sets up a pressure difference near the conductor. By using it, you can rotate the dome.
- ✓ Attach two identical sharp or rounded conductors tangentially (not radially) to the collector dome at its equator (along the seam) on opposite sides and in opposite directions. Conductors can be attached with clay or tape. Observe how pressure differences in the vicinity of these conductors exert torques on the collector dome which begin to rotate slowly but steadily. The dome continues to rotate as long as the generator is running. The mass of the dome is in excess of 200 grams excluding masses of the conductors.
- ✓ The fact, therefore, that the dome will rotate solely due to the electric wind that is generated is a testimony to the strength of that electric wind.
- **Carry the Electric Wind:** In this experiment you bring the wind to the candle instead of bringing the candle to the generator to observe its effect on the flame. Prepare a large darning needle by securely attaching a well-insulated copper wire in the needle's eye. Attach the other end of the wire to the collector dome with transparent tape. Carry the needle as far as the wire will allow you to carry it. Place it near a candle and watch the electric wind (emanating from the needle's tip) deflect the flame or turn a vane or rotate a spinner.

Demonstration:(3) Lightning


- Lightning, an awesome natural phenomenon, is an electrical discharge between clouds and the ground.
- Create it in miniature with a Van de Graaff Generator due to the buildup of positive electrical charges on the dome.
- Bring a rounded object (metallic, for best results, such as a mixing bowl or juice can or **10-074 Discharge Wand**) near the dome.
- You may wish to wear a glove or use a dry towel to hold the objects as you approach the dome to minimize the likelihood of receiving a shock.
- The discharge that occurs between rounded object and collector dome is accompanied by a crackling sound and can be made brighter and more frequent by bringing the rounded object closer (from 2" to $\frac{1}{2}$ " away.)
- If you withdraw the rounded object, the discharges become feeble and less frequent and may be seen only in a darkened room.
- You should hear intermittent crackling sounds and see feeble sparks in darkness.

(4) St. Elmo's Fire:

Electrical discharges from clouds to the earth are of three different types.





- **Point Discharge:** There is no visible light or sound. Point discharges are responsible for the bulk of discharge between clouds and ground.
- **Corona Discharge:** It is accompanied by visible light but no audible sound. This is known as St. Elmo's Fire.
- **Lightning Discharge:** This is accompanied by blinding light and deafening sound.
- ✓ You can create St. Elmo's Fire in a dark room by installing a sewing needle perpendicular to the dome using a drinking straw or small plastic strip. Tape the needle to one end of the straw, hold the straw by the other end and press it lightly against the dome. (The object, of course, is not to obtain a shock as you approach the dome with your hand.)
- ✓ A small but significant glow or "fire" appears at the tip of the needle. St. Elmo's Fire can also be created by attaching a 3' long electrical wire (not solid, but stranded) to the eye of a sewing needle.
- ✓ As the strands are passed across the eye, fold and twist them with pliers to join the needle solidly to the wire's end. Connect the other end of this wire to the ground connector on the base of your Van de Graaff. (This procedure will not work if your receptacle has only two flat holes.)
- ✓ Now tie the needle perpendicularly to one end of a drinking straw using cord or tape. Hold the far end of the straw and bring the needle close to the dome to watch the "fire" glow. With this method, you can study the effect of distance on the glow. The glow will be stronger in the vicinity of the dome. As distance increases, the glow dims. Determine the "firing distance" - the distance over which the glow is visible.

(5) Lighting

- You can light a variety of light emitting devices with your Van de Graaff Generator, for example – incandescent (filament type) light bulbs, fluorescent tubes or lamps, gas filled tubes, old radio tubes, even tiny neon tubes. For best results, do these experiments in a darkened room or at night.
- Bring your bulb toward the dome as the generator is operating. You may wish to make a nonconducting holder for the light bulb to avoid receiving a shock as you approach the dome. The outside glass surface nearest the dome acquires negative charge by induction. The charge builds up on the glass surface to discharge intensity. As discharge occurs, negative charges rush through the entire bulb, lighting it up for the duration of the discharge.
- Experiment with distances between bulb and dome. The bulb will light even when 2' away from the dome. Here, discharges will be stronger but the intervals between them will be longer. The light bulb will also glow more brightly. When you bring the bulb nearer, the discharges are more frequent but the light is dimmer. The bulb touches the dome, the light may be continuous (or flickering) but the intensity is low.
- Household (incandescent) bulbs will glow with purple light. Other gas-filled tubes will glow with the characteristic lights of the respective gases.

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
[Levitron-2](#)


<http://www.levitron.com/physics.html#intro>

- [How does the LEVITRON work](#)
- [What holds the top up?](#)
- [Why does it need to spin?](#)
- [Why doesn't the top slip sideways?](#)
- [Why is the weight so critical?](#)
- [Why does the top eventually fall?](#)
- [Is the LEVITRON Principle used elsewhere?](#)

How does the Levitron work?

What hold the top up?

- The 'antigravity' force that repels the top from the base is magnetism.
 - Both the top and the heavy slab inside the base box are magnetized, but oppositely.
 - Think of the base magnet with its north pole pointing up, and the top as a magnet with its north pole pointing down (fig 1).
 - The principle is that two similar poles (e.g., two norths) repel and that two opposite poles attract, with forces that are stronger when the poles are closer.
 - There are four magnetic forces on the top: on its north pole, repulsion from the base's north and attraction from the base's south, and on its south pole, attraction from the base's north and repulsion from the base's south.
 - Because of the way the forces depend on distance, the north-north repulsion dominates, and the top is magnetically repelled.
 - It hangs where this upward repulsion balances the downward force of gravity, that is, at the point of equilibrium where the total force is zero.
- Dr. Michael V. Berry of the University of Bristol. Dr. Berry, inspired by this recognition, published a thorough exposition of the physics of the LEVITRON's operation.

The LEVITRON[®]

A Challenge To Master; A Miracle To Behold!

Read carefully and you too will find the PATH TO SUCCESS.

The key elements you need to master in order to successfully levitate the top are detailed in the following sections:

- Section I How to spin the magnetic top on the lifter plate over the center of the magnetic base.
- Section II How to determine the correct amount of weight that the top requires to float.
- Section III How to orient the magnetic field of the base magnet.

I. THE SPIN "First let's learn to spin the magnetic top over the center of the very strong base magnet."

Before trying to spin the top, take off all of the weights (washers) except for the large brass washer. Secure the washer to the spindle with one of the rubber grommets provided (black "O-Ring").

Next try to spin the top using a slight downward pressure on the lifter plate, keeping in mind that the top wants to flip over until its spinning at least 18 RPS (Revolutions Per Second). For this reason it is important that during the "spin-up" the spindle be grasped between the thumb & forefinger tightly (pinched) and as low as possible on the spindle. The tip of your fingers should actually be touching the brass washer and parallel to the shaft of the spindle. **NOTE!** Keep your elbow raised! By keeping your elbow raised, you will also find it much easier to keep the shaft of the spindle straight up and down during the spin.

Ideally you would like to have the top spin smoothly with the spindle perfectly upright but a spin in which the spindle wobbles can also be a successful spin! By raising the lifter plate about 1/4" immediately after the "spin-up" even a wobbly spin may straighten out and become a successful spin.

If, after reading these instructions you are still not able to master the spin, you can contact your dealer and ask about the Levitron Self Starter. This is a low cost battery operated device which when placed over the tip of the spindle can be used to "spin-up" the top.

II. THE WEIGHT "Now let's learn how to determine how much weight the top will require to levitate."

The top comes with the following assortment of weight-adjustment washers which are used to balance the downward gravitational force against the upward magnetic force:

2	Large Brass	3 grams each
3	Small Brass	1 gram each
2	Red Plastic	0.4 gram each
1	Green Plastic	0.2 gram each
2	Yellow Plastic	0.1 gram each
2	Black Rubber Grommets ("O-Rings")	used to secure the weights to the spindle

The top gains stability by adding weight. The goal is to make the top float as stable as possible. This is accomplished by adding as much weight as possible to the top and still have it levitate. Of course if too much weight is added, the top will not float.

The top gains stability by adding weight. The goal is to make the top float as stable as possible. This is accomplished by adding as much weight as possible to the top and still have it levitate. Of course if too much weight is added, the top will not float.

Now, let's try to float the top. **Start with one large brass washer (3 grams) secured on the spindle:**

Step 1. Spin the top on the lifter plate over the center of the base magnet. (Practiced in Section I.)

Step 2. Raise the lifter plate approximately 1" above the surface of the base magnet. At this point continue lifting but very slowly (appx. 1/4" per 5 second interval). One customer commented that this lifting process is similar to carefully letting a car's clutch out! If the top is correctly weighted and the base magnet is properly oriented (Section III), the top will gently lift itself off the lifter plate and float in space when the top is raised to approximately 1.25 inches above the surface of the base magnet (or appx. 2.5 inches for the Super Levitron). However, chances of getting the weight and orientation correct on the first try is rare.

Step 3. Learn how to interpret the results of your levitation attempt. First ask yourself, "was the top light enough to lift off the plate as the critical height was approached or did the top just slide off the side of the plate as I continued lifting?" If the top slid off the side of the plate, then remove some of the weight from the top (for example you might remove the 3 gram weight & replace it with 2-1 gram weights) and go back to step one. If the top "jumped upward" off the plate it is too light; add more weight & try again.

By using this trial & error method, eventually you will determine the correct weight within 0.1 gram, i.e. 1 small yellow washer. The top will float only if its weight lies within a very narrow (± 0.1 gram) range of a "correct weight". You will find this correct weight by successive trials as described above.

III. THE ORIENTATION OF THE BASE **"Finally! Let's orient the magnetic field lines of the base magnet!"**

The final step in successful levitation requires that the magnetic field produced by the base magnet be precisely vertical. To achieve this it is sometimes necessary to insert the shims provided under one or even two edges of the base. Because we cannot see these magnetic field lines we must determine any adjustments by observing how the top reacts to the base.

As you experiment trying to find the correct weight for the top, notice whether the top tends to fall repeatedly to the same side. For example, if each time the top jumps off the lifter plate, it falls toward you and to your right, then the front right sides of the base need to be slightly raised. Two wedges (shims) are provided for this purpose. Slide the wedge slightly under the base at the midway point of the side or sides being raised. Then spin the top and lift again. Each time the procedure is repeated you must watch the top closely as it leaves the lifter plate to help determine what further adjustments of the wedges are needed.

REMEMBER! No matter how well you orient the base magnet, the top will not float stable unless it is carrying the correct amount of weight. (Refer to Step II.) **SO!** It is important to have the top at or near the correct floating weight **FIRST** before attempting final adjustment of the base magnet using the wedges.

FINAL NOTE . . . Once the correct weight for stable levitation has been found, only minor adjustments to that weight will be required for levitation in the future. Most of the weight variations are the result of minor changes in the strength of the magnetic material due to temperature changes. (Magnets tend to lose strength as they heat up.)

MAY THE FORCE OF LEVITATION BE WITH YOU AS YOU ATTEMPT TO MASTER the LEVITRON! GOOD LUCK!

Why does it need to spin?

- To prevent the top from overturning.
- As well as providing a force on the top as a whole, the magnetic field of the base gives a torque tending to turn its axis of spin.
- If the top were not spinning, this magnetic torque would turn it over.
- Then its south pole would point down and the force from the base would be attractive - that is, in the same direction as gravity - and the top would fall.
- When the top is spinning, the torque acts gyroscopically and the axis does not overturn but rotates about the (nearly vertical) direction of the magnetic field.
- This rotation is called **precession** (fig 2).
- With the LEVITRON the axis is nearly vertical and the precession is visible as a shivering that gets more pronounced as the top slows down.
- The effectiveness of spin in stabilizing a magnetically supported top such as the LEVITRON was discovered by Roy M. Harrigan (4).

Why doesn't the top slip sideways?

- For the top to remain suspended, equilibrium alone is not enough.
 - ⇒ The equilibrium must also be stable, so that a slight horizontal or vertical displacement produces a force pushing the top back toward the equilibrium point.
 - ⇒ For the LEVITRON stability is difficult to achieve.
- It depends on the fact that as the top moves sideways, away from the axis of the base magnet, the magnetic field of the base, about which the top's axis precesses, deviates slightly from the vertical (fig. 2). If the top precessed about the exact vertical, the physics of magnetic fields would make the equilibrium unstable.
- Because the field is so close to vertical, the equilibrium is stable only in a small range of heights - between about 1.25 inches and 1.75 inches above the center of the base. (between 2.5 and 3.0 inches for Fascinations' new Super LEVITRON).
- The Earnshaw theorem is not violated by the behavior of the LEVITRON. That theorem states that no static arrangements of magnetic (or electric) charges can be stable, alone or under gravity. It does not apply to the LEVITRON because the magnet (in the top) is spinning and so responds dynamically to the field from the base.

Why is the weight critical (and why must it be adjusted so often)?

- The weight of the top and the strength of magnetization of the base and the top determine the equilibrium height where magnetism balances gravity.
- This height must lie in the stable range. Slight changes of temperature alter the magnetization of the base and the top. (as the temperature increases, the directions of the atomic magnets randomize and the field weakens).
- Unless the weight is adjusted to compensate, the equilibrium will move outside the stable range and the top will fall. Because the stable range is so small, this adjustment is delicate - the lightest washer is only about 0.3% of the weight of the top.

Why does the top eventually fall?

- The top spins stable in the range from about 20 to 35 revolutions per second (rps). It is completely unstable above 35-40 rps and below 18 rps.
- After the top is spun and levitated, it slows down because of air resistance. After a few minutes it reaches the lower stability limit (18 rps) and falls.
- The spin lifetime of the LEVITRON can be extended by placing it in a vacuum. In a few vacuum experiments that have been done the top fell after about 30 minutes.
- Why it does so is not clear; perhaps the temperature changes, pushing the equilibrium out of the stable range; perhaps there is some tiny residual long-term instability because the top is not spinning fast enough; or perhaps vibrations of the vacuum equipment jog the field and gradually drive the precession axis away from the field direction.
- Levitation can be greatly prolonged by blowing air against an appropriately serrated air collar placed around the top's periphery so as to maintain the spin frequency in the stable range.
- Recently a LEVITRON top was kept rotating for several days in this way. But the most successful means to prolong the top's levitation is with Fascinations' new PERPETUATOR ? an electro-magnetic pulsed device which can keep the top levitating for many days or even weeks.


Is the LEVITRON Principle used elsewhere?

- In recent decades, microscopic particles have been studied by trapping them with magnetic and/or electric fields. There are several sorts of traps.
- For example, neutrons can be held in a magnetic field generated by a system of coils. Neutrons are spinning magnetic particles, so the analogy of such a neutron trap with the LEVITRON is close.

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
Ghostly Optical Illusion

- Light to Create a 3D Optical Illusion.
- You will be dazzled and amazed when your hand passes right through the floating object hovering above the Magic Mirage - Optical Illusion Generator.
- Any object placed inside the Mirage is reproduced in strikingly precise detail and floats in 3D space right before your eyes.
- This phantom reality is produced by **two precision-polished parabolic mirrors with improved Microglass Mirror Protection**.
- Comes complete and ready to use with instructions and scientific explanation.
- A great source of fun and mystery.


This phantom reality is produced by two precision-polished parabolic mirrors with improved Microglass Mirror Protection. Comes complete and ready to use with instructions and scientific explanation.

A great source of fun and mystery.


Features:

- Acrylic 9" in diameter
- No Electricity
- Amazingly simple design yet ahh-inspiring!
- Optic Mirage

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
Only electric field component along polarizing direction of polarizing sheet is passed (transmitted), the perpendicular component is blocked (absorbed).


Polarized Light

(a)

Unpolarized or randomly polarized light has its instantaneous polarization direction vary randomly with time.

One can produce unpolarized light by the addition (superposition) of two perpendicularly polarized waves with randomly varying amplitudes.


(b)

If the two perpendicularly polarized waves have fixed amplitudes and phases, one can produce different polarizations such as circularly or elliptically polarized light.

Intensity of Transmitted Polarized Light

One-half Rule

Intensity of transmitted light, unpolarized incident light: $I = \frac{1}{2} I_0$

Since only the component of the incident electric field E parallel to the polarizing axis is transmitted

$$E_{\text{transmitted}} = E_y = E \cos \theta$$

Intensity of transmitted light, polarized incident light: $I = I_0 \cos^2 \theta$

Cosine-sqaure rule: For unpolarized light, θ varies randomly in time

$$I = \left(I_0 \cos^2 \theta \right)_{\text{avg}} = I_0 \left(\cos^2 \theta \right)_{\text{avg}} = I_0 \frac{1}{2}$$

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
Maxwell's Rainbow

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The wavelength/frequency **range** in which **EM waves** (light) are visible is only a tiny fraction of the entire electromagnetic **spectrum**

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
Euler's disk

- The disk is a chrome-plated steel disk with nine pieces of magnetized holographic foil.
- A concave mirror base is the setting for the disk to spin like a coin—except it seems like it will never stop!!
- Just give it a spin and gravity does the rest!!

Operation:

- Observing the images and sound as the disk spins faster and louder.
- Want an increased visual effect? Shine a flashlight on the disk while it is spinning in a darkened room!

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Hydrogen Fuel Car- H Racer

The smallest hydrogen fuel cell car in the world!

Solar panel

Hydrogen Storage Station

Racing Car

Hydrogen Fuel Cell

The car uses a real fuel cell and its own on-board hydrogen storage system. It does not need batteries! The Hydrogen Station will provide the H-racer with an unlimited supply of clean energy.

To create free hydrogen fuel at the flick of a switch, just add water to the station's tank! Fueling is animated by a special blue light display.

- The **H-racer** is a futuristic toy car that contains one of the most exciting and advanced technologies of the 21st century. This car operates on 100% clean fuel produced by a miniature solar-powered hydrogen refueling station that converts water to hydrogen using energy captured from the sun.
- Cars running on clean and renewable fuel are the dream of many of today's world leaders, engineers and scientists seeking to eliminate mankind's reliance on fossil fuels in favor of climate friendly energy resources.

- With new advances in technology, hydrogen is on its way to becoming the world's next fuel. Hydrogen offers many important advantages: it is non-toxic, renewable, clean to use, and the most abundant element in our universe. And by using fuel cell technology to convert hydrogen to electricity without any combustion, the technology is a significant solution to many of our global energy and environmental problem.
- No combustion occurs inside a fuel cell. The only exhaust resulting from hydrogen fuel cell cars is pure water. Fuel cell cars that use hydrogen as a fuel are also known as “zero emission vehicles.”
- Today, many of the world's automotive companies including Toyota, GM, Ford, Honda, and Daimler-Chrysler, are developing hydrogen fuel cell vehicles with the hope of introducing this technology to the public in the near future.
- The **H-racer** is the working miniature version of what is being developed in real-size cars of the future.
- This palm-size fuel cell car contains an onboard hydrogen storage tank, a fuel cell system connected to the car's electric motor, and a hydrogen refueling system linking the car's storage tank to an external hydrogen refueling station.
- Given its small size, the **H-racer** is also very safe as only tiny quantities of hydrogen are sufficient to power the car.